

Eisner™

A division of Multigate  
Medical Products Pty Ltd



## Forcep Inserts

Eisner™ laparoscopic  
instruments deliver  
unmatched durability,  
comfort and efficiency  
to support a seamless  
performance every time.

1993

POWERING  
LAPAROSCOPIC  
SURGERY  
SINCE 1993



**DESIGN**

Our instruments are designed for optimal balance, feel, comfort and responsiveness to deliver a seamless connection between the hand and the surgical site.

**SAFETY & RELIABILITY**

The craftsmanship and the quality of materials we select combine to deliver high performance, electro-safe instruments that deliver unrivalled durability, longevity and reliability.



- **Easy to re-process** Our instruments are designed to simplify re-processing.

- **Clinical performance** Our hand-held instruments are modular allowing you to build cost effective instrument trays whilst maximising patient safety and clinical performance.





## INSERTS - FORCEPS

Eisner™ inserts are tough performers designed to deliver high fidelity control without compromising on their durability, lifespan and reliability.

- **Easy to manage** Our inserts are designed to self-open for sterilisation.

- **Durable** No weld design ensures long-term performance.



**Working Length** Inserts available in standard and bariatric lengths.

**Performance** Broad range of precision jaws to support all types of procedure.

## Inserts - Forceps (5mm)

We offer a broad range of jaw types in both standard and bariatric lengths to cover all minimally invasive procedures.

All our inserts are fully compatible and inter-changeable and designed for durable, reliable performance.

Biopsy	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>5mm Biopsy Forceps</b>	400082	401082	11mm Jaw	Double Action Sharp Oval Cups	1
	<b>5mm Biopsy Pin Forceps</b>	400521	401521	10mm Jaw	Single Action Sharp Oval Cups 2 Pins	1
	<b>5mm Biopsy Punch Insert</b>	400100	401100	11mm Jaw	Single Action Thru-Cutting Mechanism	1
	<b>5mm Blakesley Biopsy Forceps</b>	400086	401086	10mm Jaw	Single Action Sharp Elongated Fenestrated Oval Cups	1
	<b>5mm Minibite Biopsy Forceps</b>	400080	401080	8mm Jaw	Double Action Sharp Short Oval Cups 2 Pins	1

Ovary/Tubal	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>5mm Blunt Nose Tubal Forceps</b>	400046	401046	20mm Jaw	Double Action Ovary/Tubal Forceps Non Serrated Inner Jaw	1
	<b>5mm Tubal Forceps</b>	400045	401045	16mm Jaw	Double Action Fine Horizontal Serrations	1

## INSERTS - FORCEPS (5MM)

Clamps	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>5mm Debakey Forceps</b>	400102	401102	20mm Jaw	Double Action 2 into 3 Debakey Serrations Slim Line Profile	1
	<b>5mm Debakey Forceps</b>	400103	401103	30mm Jaw	Double Action 2 into 3 Debakey Serrations Slim Line Profile	1
	<b>5mm Debakey Forceps</b>	400104	401104	40mm Jaw	Double Action 2 into 3 Debakey Serrations Slim Line Profile	1
	<b>5mm Dorsey Fenestrated Forceps</b>	400287	401287	45mm Jaw	Double Action Heavy Horizontal Serrations Long Single Fenestration	1
	<b>5mm Hunter Jaeger Forceps</b>	400012	401012	21mm Jaw	Double Action Debakey Serrations Long Single Fenestration	1
	<b>5mm Marques DA Double Fen. Forceps</b>	400311	401311	25mm Jaw	Double Action Medium Horizontal Serrations Double Fenestration	1
	<b>5mm Marques DA Triple Fen. Forceps</b>	400313	401313	35mm Jaw	Double Action Medium Horizontal Serrations Triple Fenestration	1
	<b>5mm Marques SA Double Fen. Forceps</b>	400310	401310	25mm Jaw	Single Action Medium Horizontal Serrations Double Fenestration	1

## INSERTS - FORCEPS (5MM)

Clamps (continued)	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>5mm Marques SA Triple Fen. Forceps</b>	400312	401312	35mm Jaw	Single Action Medium Horizontal Serrations Triple Fenestration	1
	<b>5mm Marques Upcurve Forceps</b>	400314	401314	25mm Jaw	Single Action Medium Horizontal Serrations Double Fenestration	1

Curved Vascular Clamps	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>5mm Crafoord 35 Clamp Forceps</b>	400204	401204	35mm Jaw	Single Action Atraumatic 1 into 2 Debakey Serrations	1
	<b>5mm Crafoord 50 Clamp Forceps</b>	400206	401206	50mm Jaw	Single Action Atraumatic 1 into 2 Debakey Serrations	1
	<b>5mm Debakey Curved 50 Clamp Forceps</b>	400220	401220	50mm Jaw	Single Action Atraumatic 1 into 2 Debakey Serrations	1
	<b>5mm Debakey S-shape 50 Clamp Forceps</b>	400221	401221	50mm Jaw	Single Action Atraumatic 1 into 2 Debakey Serrations	1
	<b>5mm Satinsky 40 Clamp Forceps</b>	400202	401202	40mm Jaw	Single Action Atraumatic 1 into 2 Debakey Serrations	1
	<b>5mm Satinsky 50 Clamp Forceps</b>	400200	401200	50mm Jaw	Single Action Atraumatic 1 into 2 Debakey Serrations	1

## INSERTS - FORCEPS (5MM)

Dissectors	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>5mm Dolphin Cup Forceps</b>	400015	401015	15mm Jaw	Double Action Medium Horizontal Serrations with Cup	1
	<b>5mm Dolphin Forceps</b>	400021	401021	14mm Jaw	Double Action Medium Horizontal Serrations	1
	<b>5mm Duckbill Rounded Forceps</b>	400087	401087	10mm Jaw	Double Action Parallel/Aligned Jaw Closure Medium Horizontal Serrations	1
	<b>5mm Geegee Forceps</b>	400030	401030	38mm Jaw	Single Action Medium Horizontal Serrations	1
	<b>5mm Kelly Baby Criss Cross Forceps</b>	400733	401733	15mm Jaw	Double Action Criss Cross Serrations	1
	<b>5mm Kelly Criss Cross Forceps</b>	400732	401732	20mm Jaw	Double Action Criss Cross Serrations	1
	<b>5mm Maryland Debakey Forceps</b>	551712	—	29mm Jaw	Double Action 1 into 2 Debakey Serrations Fine Teeth at Tip	1
	<b>5mm Maryland Long Forceps</b>	400660	401660	24mm Jaw	Double Action Criss Cross Serrations	1

## INSERTS - FORCEPS (5MM)

Dissectors (continued)	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>5mm Maryland Standard Forceps</b>	400731	401731	20mm Jaw	Double Action Fine Horizontal Serrations	1
	<b>5mm Maryland Grasper Dissector Forceps</b>	400661	401661	17mm Jaw	Double Action Deep Horizontal Serrations	1
	<b>5mm Mixer Forceps</b>	400735	401735	13mm Jaw	Double Action Parallel/Aligned Jaw Closure Fine Criss Cross Serrations	1
	<b>5mm Pointed Forceps</b>	400016	401016	22mm Jaw	Double Action Criss Cross Serrations	1
	<b>5mm Reddick Olsen Forceps</b>	400342	401342	13mm Jaw	Double Action Medium Horizontal Serrations	1

Fenestrated Dissectors	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>5mm Bullnose Fenestrated Forceps</b>	400040	401040	10mm Jaw	Double Action Medium Horizontal Serrations Single Fenestration	1
	<b>5mm Fenestrated Dissector Forceps</b>	400088	401088	13mm Jaw	Double Action Fine Horizontal Serrations Single Elongated Fenestration	1

## INSERTS - FORCEPS (5MM)

Fenestrated Dissectors (continued)	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>5mm Fenestrated Strong Forceps</b>	400132	401132	8mm Jaw	Double Action Medium Horizontal Serrations Single Elongated Fenestration	1
	<b>5mm Johans Fenestrated Forceps</b>	400285	401285	17mm Jaw	Double Action Medium Horizontal Serrations Spacious Fenestration	1
	<b>5mm Johans Long DA Fen. Forceps</b>	400291	401291	40mm Jaw	Double Action Super Fine Horizontal Serrations Long Spacious Fenestration	1
	<b>5mm Johans Long SA Fen. Forceps</b>	400288	401288	40mm Jaw	Single Action Super Fine Horizontal Serrations Long Spacious Fenestration	1
	<b>5mm Johans Loop DA Forceps</b>	400289	401289	21mm Jaw	Double Action Medium Horizontal Serrations Spacious Fenestration	1
	<b>5mm Johans Loop SA Forceps</b>	400284	401284	20mm Jaw	Single Action Medium Horizontal Serrations Spacious Fenestration	1
	<b>5mm Wave Fenestrated Forceps</b>	400052	401052	14mm Jaw	Double Action Fine Horizontal Serrations Single Long Fenestration	1

## INSERTS - FORCEPS (5MM)

Extractors	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>5mm Cats Paw Tenaculum Forceps</b>	400048	401048	20mm Jaw	Double Action 2 into 3 End Teeth 2 Pair Outer Side Hook Teeth	1
	<b>5mm Davis Claw (Serrated Tooth) Forceps</b>	400717	401717	22mm Jaw	Double Action 2 Into 3 End Teeth Deep Horizontal Grooves	1
	<b>5mm Eistige Forceps</b>	400013	401013	22mm Jaw	Double Action Deep Cup Rippled Edges	1
	<b>5mm Eistige SA Forceps</b>	400033	401033	22mm Jaw	Single Action Deep Cup Rippled Edges	1
	<b>5mm Gripper Forceps</b>	400017	401017	21mm Jaw	Double Action Coarse Serrations Long Fenestration	1
	<b>5mm Interlocker Forceps</b>	400217	401217	22mm Jaw	Double Action Interlocking Smooth Raised Teeth Long Fenestration	1
	<b>5mm Ripple Claw (St Andrews) Forceps</b>	400715	401715	23mm Jaw	Double Action 2 into 4 Aggressive End Teeth	1
	<b>5mm Tenaculum Forceps</b>	400364	401364	22mm Jaw	Double Action Super Fine Claw Tip Smooth Inner Jaw	1

## INSERTS - FORCEPS (5MM)

Graspers	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>5mm Alligator Retrograde Forceps</b>	400010	401010	15mm Jaw	Double Action Strong Retrograde Serrations	1
	<b>5mm Allis 35 Forceps</b>	400197	401197	35mm Jaw	Double Action 4mm Wide Double Row Platform Smooth Inner Jaw	1
	<b>5mm Allis Forceps</b>	400089	401089	13mm Jaw	Double Action Serrated End Platform Smooth Inner Jaw	1
	<b>5mm Allis SA Forceps</b>	400011	401011	13mm Jaw	Single Action Serrated End Platform Smooth Inner Jaw	1
	<b>5mm Allis Long Forceps</b>	400330	401330	28mm Jaw	Double Action Fine Toothed Tip Smooth Inner Jaw	1
	<b>5mm Babcock Duval Forceps</b>	400130	401130	25mm Jaw	Double Action Horizontal Serrations on Medium Platform Elongated Fenestration	1
	<b>5mm Babcock SA Forceps</b>	400251	401251	12mm Jaw	Single Action Fine Horizontal Serrations Long Fenestration	1
	<b>5mm Babcock Super Fine Paddle Forceps</b>	400250	401250	30mm Jaw	Double Action 5mm Platform with Horizontal Serrations Smooth Inner Jaw	1

## INSERTS - FORCEPS (5MM)

Graspers (continued)	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>5mm Babcock Traditional Forceps</b>	400252	401252	20mm Jaw	Double Action Distal End Fenestration Small Platform with Horizontal Serrations	1
	<b>5mm Bonnies Forceps</b>	400009	401009	15mm Jaw	Double Action 1 into 2 Teeth Fine Horizontal Serrations	1
	<b>5mm Bullnose Flat Face Forceps</b>	400355	401355	13mm Jaw	Double Action Medium Horizontal Serrations	1
	<b>5mm Bullnose Forceps</b>	400127	401127	13mm Jaw	Double Action Medium Horizontal Serrations	1
	<b>5mm Bullnose Long Forceps</b>	400002	401002	22mm Jaw	Double Action Medium Horizontal Serrations	1
	<b>5mm Bullnose Short Jaw Forceps</b>	400001	401001	11mm Jaw	Double Action Medium Horizontal Serrations	1
	<b>5mm Crocodile Forceps</b>	400019	401019	17mm Jaw	Double Action Precision Pyramidal Teeth Small Triple Fenestration	1
	<b>5mm Kocher Forceps</b>	400004	401004	15mm Jaw	Double Action 2 into 3 End Teeth Fine Horizontal Serrations	1

## INSERTS - FORCEPS (5MM)

Graspers (continued)	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>5mm McKearnan Medium Forceps</b>	400136	401136	15mm Jaw	Double Action Traditional McKearnan Style	1
	<b>5mm McKearnan Standard Forceps</b>	400025	401025	20mm Jaw	Double Action Traditional McKearnan Style	1
	<b>5mm McKearnan Long Forceps</b>	400137	401137	22mm Jaw	Double Action McKearnan Style Long Single Fenestration	1
	<b>5mm Rat Tooth Forceps</b>	400036	401036	12mm Jaw	Double Action 2 into 4 Teeth Medium Horizontal Serrations	1
	<b>5mm Pin Grasper Forceps</b>	400042	401042	12mm Jaw	Double Action 1 into 2 End Teeth Smooth Inner Jaw	1
	<b>5mm Tungsten Carbide (TC) Grasper Forceps</b>	400029	401029	16mm Jaw	Double Action Tungsten Carbide Jaw Pyramidal Serrations	1
	<b>5mm Cobra Forceps</b>	400201	401201	13mm Jaw	Double Action 2 into 4 Sharp End Teeth	1
Grasper Scissors	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>5mm Knot Pusher Scissor Insert</b>	400331	401331	16mm Jaw	Double Action Knot Pusher Tip and Scissors in Unique Molding	1

## Inserts - Forceps (10mm)

The density of our 10mm insulation tubes is tripled to minimize splitting, cracking and deformity and their single piece design with no welded joints delivers heavy duty, reliable performance every time.

Biopsy	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>10mm Spoon Forceps</b>	400720	401720	27mm Jaw	Single Action Sharp Large Cups Rounded End Profile	1
Clamps	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>10mm Nussbaum Intestinal Forceps</b>	400781	401781	90mm Jaw	Double Action Full Length Vertical Serrations Additional Fine Cross-Serrations on Leading Section	1
Dissectors	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>10mm Geegee Forceps</b>	400076	401076	52mm Jaw	Single Action Medium Horizontal Serrations	1
	<b>10mm Mixer 90° Forceps</b>	400739	401739	28mm Jaw	Double Action Medium Horizontal Serrations	1
	<b>10mm Overholt Criss Cross Forceps</b>	400741	401741	25mm Jaw	Double Action Medium Criss Cross Serrations	1
	<b>10mm Overholt-Mixer Stepdown Forceps</b>	400737	401737	30mm Jaw	Double Action Deep Horizontal Serrations	1
	<b>10mm Right Angle Forceps</b>	400736	401736	24mm Jaw	Double Action Fine Criss Cross Serrations	1

## INSERTS - FORCEPS (10MM)

Extractors	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>10mm Claw (Mother in-law) Forceps</b>	400710	401710	30mm Jaw	Single Action Aggressive 2 into 3 End Teeth Smooth Inner Jaw	1
	<b>10mm Tenaculum Forceps</b>	400753	401753	35mm Jaw	Double Action Fine Claw Tip Smooth Inner Jaw	1
Graspers	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	<b>10mm Anvil Grasping Forceps</b>	400295	401295	30mm Jaw	Double Action Specially Shaped to Hold the Anvil of a Circular Stapler	1
	<b>10mm Babcock Forceps</b>	400078	401078	30mm Jaw	Double Action Medium Platform with Fine Vertical Serrations V-shaped Fenestration	1
	<b>10mm Collin-Duval Grasping Forceps</b>	400084	401084	32mm Jaw	Double Action Criss Cross Serrations on Platform Long Broad Fenestration	1
	<b>10mm Duval Long Forceps</b>	400280	401280	50mm Jaw	Double Action Deep Vertical Serrations on Platform Short V-shaped Fenestration	1
	<b>10mm Pennington-Collin Forceps</b>	400085	401085	38mm Jaw	Double Action Deep Horizontal Serrations Spacious V-shaped Fenestration	1

## SUPPORT

With Eisner™ you don't get a call centre, just a dedicated team of specialists ready to tailor your Eisner™ experience.

### Experience the difference

#### Book an obligation free trial

Our team can organise an obligation free, trial assessment of your instruments so you can see and feel the difference before you buy.

### Trade-in / replacement program

#### Manage your tray upgrades

If you are considering decommissioning an instrument our trade-in/replacement service will assist you to address one off replacements or the gradual upgrade of instrument sets. We accept Eisner™ and non-Eisner™ instruments.

### Care

#### Optimise Your Tray Performance

When you choose Eisner™ products our team are available to provide free comprehensive in-service training, support and tray assessments to guide servicing and repair.





# Enhance your every move

**Warranty** | Eisner™ offers a 12 month warranty against manufacturing defect on all metal components of the equipment. Damage and / or breakages caused by improper use or abuse of the equipment are not covered by warranty. Normal wear and tear is not covered. Warranty is null and void where components are mixed with non-Eisner™ instruments.

**Service & Repairs** | Instruments returned to us for repair/replacement will be assessed and a written quotation provided for you to make a decision on further action. In most cases for instruments beyond repair we can offer a trade-in replacement. Eisner™ keeps factory original spare parts in stock. All repairs should be returned to Eisner™ for repair in order to preserve the integrity of your instruments.

**Recommended Lifespan of Instruments** | Whilst Eisner™ laparoscopic instruments will function for many years, we recommend replacement prior to any obvious signs of wear and stress which could lead to an issue during surgery. We recommend a lifespan of 3 years on our handles and inserts and suggest that after that time, they are traded in on a new component. Insulated Tubes: By virtue of the fact that insulated tubes will be replaced once the insulation is damaged or compromised in any way, it is not envisioned that the lifespan of a sheath will be beyond 12 months. Instruments should be taken out of circulation immediately if it is found that identifying labelling is no longer visible.

**Eisner**<sup>™</sup>

A division of Multigate  
Medical Products Pty Ltd

Pacifica Business Park  
Unit 18, 97 Canterbury Road  
Kilsyth VIC 3137  
Australia  
T +61 (0)3 9728 4770  
F +61 (0)3 9728 6770  
eisner@multigate.com.au  
[www.eisner.com.au](http://www.eisner.com.au)

For an obligation  
free trial or for more  
information call  
us today on:  
(03) 9728 4770  
[eisner.com.au](http://eisner.com.au)